

## Rhinorolls



<b>COMPOSITION</b>	<p>100% polyamide tubular nose dressing containing an absorbent roll of:</p> <ul style="list-style-type: none"> <li>• 25% cotton</li> <li>• 25% viscose</li> <li>• 50% cellulose</li> </ul>
<b>RHINOROLL</b>	<ul style="list-style-type: none"> <li>• High elasticity thanks to the polyamide</li> <li>• One size fits all: +/- 42 cm long, absorbent roll of 8 cm long, Ø 15mm</li> <li>• With ear-loops</li> <li>• Used after nose surgery to absorb secretion emission</li> <li>• Can be used as support for decongestant treatment during common respiratory diseases</li> <li>• Non sterile</li> <li>• Shelf life : 5 years from production date</li> </ul>

### Packaging

References	Description	Pieces/Box	Weight per box	Boxes/Export box	Export boxes/Pallet
RRA20	Rhinorolls	20	+/- 120g	100	10

### Comments :

Class I medical device, single use product  
 Packed in cardboard boxes 23x9x6 cm to protect from dust  
 Export box sizes: 59x38x72 cm  
 General tolerance: +/- 10%  
 Should be stored in a clean and dry place, out of direct sunlight

**We are ISO 13485 certified and our products are CE marked**

