

Astringent solution



| | |
|----------------------------|---|
| COMPOSITION | <ul style="list-style-type: none"> • Octadecahydrated aluminium sulphate, PH.Eur. (25%) • Distilled Water (75%) • Without epinephrine : no risk of circulatory disorders |
| ASTRINGENT SOLUTION | <ul style="list-style-type: none"> • Astringent and haemostatic solution used for treatment of gingival retraction together with our retraction cord, knitted or braided • For dental use only • Less sensitive to hydrolysis than Aluminium Chloride (more stable) • Colourless, non-staining • Shelf life : 3 years from production date |

Packaging :

| Reference | Description | Solution/Bottle | Bottles/Case | Cases/Box |
|-----------|--------------|-----------------|--------------|-----------|
| ALST20 | Glass bottle | 20 ml | 1 | 49 |

Comments:

Cannot be used in case of severe ulcerous gingivitis

Precaution: slightly acidic salt, rinse thoroughly in case of contact with eyes

Sold in glass bottle in order to keep the product stable

Weight of the bottle: 41,5 g

Case sizes: 30 x 30 x 63 mm

Export box sizes: 225 x 225 x 80 mm

General tolerance: +/- 10 %

Should be stored in a clean and dry place, out of direct sunlight and at room temperature

We are ISO 13485 certified and our products are CE marked

